




Product Datasheet

GALNT3 Antibody (orb1262319)

Description	GALNT3 Antibody	
Species/Host	Rabbit	Western blot analysis of GALNT3 Antibody...
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	FC, IHC-P, WB	
Immunogen	This GALNT3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 454-481 amino acids from the Central region of human GALNT3.	
Target	GALNT3	 Formalin-fixed and paraffin-embedded hum...
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	 Flow cytometric analysis of MCF-7 cells ...
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	73 kDa	
Uniprot ID	Q14435	
NCBI	Q14435	
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50	
Expiration Date	12 months from date of receipt.	