

## Product Datasheet

### MPP6 Antibody (orb1262322)

## Description

MPP6 Antibody

### Species/Host

Rabbit

### Reactivity

Human, Mouse

### Conjugation

Unconjugated

### Tested

FC, WB

### Applications

### Immunogen

This MPP6 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 360-388 amino acids from the Central region of human MPP6.

### Target

MPP6

### Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

### Note

For research use only

### Application notes

For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50

### Isotype

Rabbit Ig

### Clonality

Polyclonal

### MW

61 kDa

### Uniprot ID

[Q9NZW5](#)

### NCBI


[Q9NZW5](#)

### Dilution Range


For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50

### Expiration Date

12 months from date of receipt.



Western blot analysis of MPP6 Antibody i...



Flow cytometric analysis of 293 cells (r...