

Product Datasheet

PEX2 Antibody (orb1262334)

Description

PEX2 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications
Immunogen

This PXMP3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 165-194 amino acids from the Central region of human PXMP3.

Target

PEX2

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

35 kDa


Uniprot ID
P28328
NCBI
P28328
Dilution Range

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50


Expiration Date

12 months from date of receipt.


m.cerebellum



Western blot analysis of PXMP3 Antibody ...



Formalin-fixed and paraffin-embedded hum...



Flow cytometric analysis of WiDr cells (...)