

Product Datasheet

HOXA9 Antibody (orb1262367)

Description

HOXA9 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, WB

Applications

Immunogen

This HOXA9 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 245-272 amino acids from the C-terminal region of human HOXA9.

Target

HOXA9

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

30 kDa

Uniprot ID

[P31269](#)

NCBI

[P31269](#)


Dilution Range

For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50


Expiration Date

12 months from date of receipt.

A2058



Western blot analysis of HOXA9 Antibody ...



Flow cytometric analysis of A375 cells (...)