



Product Datasheet

MTTP Antibody (orb1262391)

Description	MTTP Antibody	
Species/Host	Rabbit	
Reactivity	Human	Western blot analysis of MTTP Antibody i...
Conjugation	Unconjugated	
Tested Applications	IHC-P, WB	Formalin-fixed and paraffin-embedded hum...
Immunogen	This MTTP antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 816-842 amino acids from the C-terminal region of human MTTP.	
Target	MTTP	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	99 kDa	
Uniprot ID	P55157	
NCBI	P55157	
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50	
Expiration Date	12 months from date of receipt.	