

Product Datasheet

TUBB1 Antibody (orb1262410)

Description
TUBB1 Antibody
Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications
Immunogen

This TBB1 antibody is generated from rabbits immunized with human TBB1 recombinant protein.

Target

TUBB1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:2000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


50 kDa

Uniprot ID
Q9H4B7
NCBI
Q9H4B7
Dilution Range


For WB starting dilution is: 1:2000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western Blot
at 1:2000
dilution Lane
1: ...



Formalin-
fixed and
paraffin-
embedded
hum...



Flow
cytometry
analysis of
MCF-7 cells
(...