



Product Datasheet

DPP10 Antibody (orb1262418)

Description	DPP10 Antibody	
Species/Host	Rabbit	 Western blot analysis of DPP10 Antibody ...
Reactivity	Human, Mouse	
Conjugation	Unconjugated	
Tested Applications	IHC-P, WB	
Immunogen	This DPP10 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 150-177 amino acids from the Central region of human DPP10.	
Target	DPP10	 DPP10 Antibody IHC analysis in formalin ...
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	91 kDa	
Uniprot ID	Q8N608	
NCBI	Q8N608	
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100	
Expiration Date	12 months from date of receipt.	