

Product Datasheet

ELAVL1 Antibody (orb1262434)


Description

ELAVL1 Antibody


Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	This ELAVL1 antibody is generated from rabbits immunized with human ELAVL1 recombinant protein.
Target	ELAVL1
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For FACS starting dilution is: 1:25For IHC-P starting dilution is: 1:25For WB starting dilution is: 1:2000
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	36 kDa
Uniprot ID	Q15717
NCBI	Q15717
Dilution Range	For FACS starting dilution is: 1:25For IHC-P starting dilution is: 1:25For WB starting dilution is: 1:2000
Expiration Date	12 months from date of receipt.



Overlay histogram showing MCF-7 cells stained with ELAVL1 antibody (green) and isotype control (blue).



Overlay histogram showing HeLa cells stained with ELAVL1 antibody (green) and isotype control (blue).



Antibody staining ELAVL1 in human testis tissue.