

Product Datasheet

COL18A1 Antibody (orb1262484)

Description

COL18A1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

FC, WB

Applications
Immunogen

This COL18A1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 237-263 amino acids from the N-terminal region of human COL18A1.

Target

COL18A1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


178 kDa

Uniprot ID
P39060
NCBI
P39060
Dilution Range


For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50

Expiration Date

12 months from date of receipt.



Western blot analysis of COL18A1 Antibod...



Flow cytometric analysis of 293 cells (r...