

Product Datasheet

PLCL1 Antibody (orb1262485)

Description

PLCL1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, WB

Applications

Immunogen

This PLCL1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 149-176 amino acids from the N-terminal region of human PLCL1.

Target

PLCL1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

123 kDa

Uniprot ID

[Q15111](#)

NCBI


[Q15111](#)

Dilution Range


For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50

Expiration Date

12 months from date of receipt.



Western blot analysis of PLCL1 Antibody ...



Flow cytometry analysis of Jurkat cells ...