


Product Datasheet

SCNN1A Antibody (orb1262503)


Description

SCNN1A Antibody


Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	This SCNN1A antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 365-391 amino acids from the Central region of human SCNN1A.
Target	SCNN1A
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:2000For FACS starting dilution is: 1:25For IHC-P starting dilution is: 1:10~50For IF starting dilution is: 1:10~50
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	76 kDa
Uniprot ID	P37088
NCBI	P37088
Dilution Range	For WB starting dilution is: 1:2000For FACS starting dilution is: 1:25For IHC-P starting dilution is: 1:10~50For IF starting dilution is: 1:10~50
Expiration Date	12 months from date of receipt.



Western Blot
at 1:2000
dilution +
Daudi ...



Overlay
histogram
showing
WiDr cells
sta...



Formalin-
fixed and
paraffin-
embedded
hum...