

Product Datasheet

IL18RAP Antibody (orb1262506)

Description

IL18RAP Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications

Immunogen

This IL18RAP antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 73-101 amino acids from the N-terminal region of human IL18RAP.

Target

IL18RAP

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

68 kDa

Uniprot ID

[O95256](#)

NCBI

[O95256](#)

Dilution Range

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50


Expiration Date

12 months from date of receipt.


250

100
75
50
37
25

Western blot analysis of IL18RAP Antibod...



Formalin-fixed and paraffin-embedded hum...



Flow cytometric analysis of k562 cells (...)