

Product Datasheet

SKIV2L2 Antibody (orb1262533)

Description

SKIV2L2 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications
Immunogen

This SKIV2L2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 237-263 amino acids from the Central region of human SKIV2L2.

Target

SKIV2L2

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

118 kDa

Uniprot ID
P42285
NCBI
P42285
Dilution Range

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Expiration Date

12 months from date of receipt.

250


130

95


72

55

Western blot analysis of SKIV2L2 Antibod...



Formalin-fixed and paraffin-embedded hum...



Flow cytometry analysis of MDA-MB231 cel...