

Product Datasheet

ACY3 Antibody (orb1262546)

Description

ACY3 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

FC, WB

Applications

Immunogen

This ACY3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 66-95 amino acids from the Central region of human ACY3.

Target

ACY3

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

35 kDa

Uniprot ID

Q96HD9

NCBI


Q96HD9

Dilution Range


For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western blot analysis of ACY3 Antibody i...



Western blot analysis of ACY3 Antibody i...



Flow cytometric analysis of HepG2 cells ...