

## Product Datasheet

### AGPAT4 Antibody (orb1262553)

## Description

AGPAT4 Antibody

### Species/Host

Rabbit

### Reactivity

Human

### Conjugation

Unconjugated

### Tested

FC, WB

### Applications

### Immunogen

This AGPAT4 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 176-202 amino acids from the Central region of human AGPAT4.

### Target

AGPAT4

### Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

### Note

For research use only

### Application notes

For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50

### Isotype

Rabbit Ig

### Clonality

Polyclonal

### MW

44 kDa

### Uniprot ID

[Q9NRZ5](#)

### NCBI

[Q9NRZ5](#)

### Dilution Range


For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50

### Expiration Date


12 months from date of receipt.

MDA-MB435

130



Western blot analysis of AGPAT4 Antibody...



Flow cytometry analysis of MDA-MB231 cel...