

Product Datasheet

ZYG11A Antibody (orb1262570)

Description

ZYG11A Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, WB

Applications

Immunogen

This ZYG11A antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 71-97 amino acids from the N-terminal region of human ZYG11A.

Target

ZYG11A

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

86 kDa

Uniprot ID

[Q6WRX3](#)

NCBI

[Q6WRX3](#)


Dilution Range

For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50


Expiration Date

12 months from date of receipt.

MCF-7
130



Western blot analysis of ZYG11A Antibody...



Flow cytometric analysis of MCF-7 cells ...