


Product Datasheet

KIR3DS1 Antibody (orb1262572)


Description

KIR3DS1 Antibody

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	This KIR3DS1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 361-387 amino acids from the C-terminal region of human KIR3DS1.
Target	KIR3DS1
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	43 kDa
Uniprot ID	Q14943
NCBI	Q14943
Dilution Range	For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50
Expiration Date	12 months from date of receipt.



Western blot analysis of KIR3DS1 Antibod...



Flow cytometric analysis of K562 cells (...)