

Product Datasheet

SCP2 Antibody (orb1262589)

Description

SCP2 Antibody

130

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications
Immunogen

This SCP2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 358-385 amino acids from the Central region of human SCP2.

Target

SCP2

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


59 kDa

Uniprot ID
P22307
NCBI
P22307
Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western blot analysis of SCP2 Antibody i...



Formalin-fixed and paraffin-embedded hum...



Flow cytometric analysis of HepG2 cells

...