



Product Datasheet

HLA-E Antibody (orb1262653)

Description	HLA-E Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	This HLA-E antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 108-135 amino acids from the Central region of human HLA-E.
Target	HLA-E
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:1000
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	40 kDa
Uniprot ID	P13747
NCBI	P13747
Dilution Range	For WB starting dilution is: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of HLA-E Antibody ...



Western blot analysis of HLA-E using rab...