

Product Datasheet

PPP3CC Antibody (orb1262665)

Description

PPP3CC Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

This PPP3CC antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 482-512 amino acids from the C-terminal region of human PPP3CC.

Target

PPP3CC

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For IHC-P starting dilution is: 1:50~100For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

58 kDa

Uniprot ID

[P48454](#)

NCBI


[P48454](#)

Dilution Range


For IHC-P starting dilution is: 1:50~100For WB starting dilution is: 1:1000

Expiration Date

12 months from date of receipt.



Formalin-fixed and paraffin-embedded hum...



Western blot analysis of anti-PPP3CC Pab...