


Product Datasheet

DUSP4 Antibody (orb1262677)


Description

DUSP4 Antibody


Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	This DUSP4 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 331-361 amino acids from the C-terminal region of human DUSP4.
Target	DUSP4
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	43 kDa
Uniprot ID	Q13115
NCBI	Q13115
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100
Expiration Date	12 months from date of receipt.



Western blot analysis of anti-DUSP4
Anti...



Western blot analysis of anti-DUSP4
Pab ...



Formalin-fixed and paraffin-embedded hum...