

Product Datasheet

PTPN3 Antibody (orb1262688)

Description

PTPN3 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

This PTPH1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 366-397 amino acids from the Central region of human PTPH1.

Target

PTPN3

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

104 kDa

Uniprot ID

[P26045](#)

NCBI

[P26045](#)


Dilution Range

For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:10~50


Expiration Date

12 months from date of receipt.

250



Antibody is used in Western blot to dete...



Formalin-fixed and paraffin-embedded hum...