

Product Datasheet

TRPM8 Antibody (orb1262698)

Description
TRPM8 Antibody
Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

This TRPM8 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 521-552 amino acids from the Central region of human TRPM8.

Target

TRPM8

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


128 kDa

Uniprot ID
[Q7Z2W7](#)
NCBI
[Q7Z2W7](#)
Dilution Range


For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western blot analysis of lysate from hum...



Western blot analysis of TRPM8 using rab...



Formalin-fixed and paraffin-embedded hum...