

Product Datasheet

BAI1 Antibody (orb1262708)

Description

BAI1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested Applications

FC, IHC-P, WB

Immunogen

This BAI1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 1537-1567 amino acids from the C-terminal region of human BAI1.

Target

BAI1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

174 kDa

Uniprot ID

[O14514](#)

NCBI


[O14514](#)

Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western blot analysis in K562 cell line ...



Western blot analysis in mouse bladder t...



BAI1 antibody immunohistochemistry analy...