



Product Datasheet

NEK3 Antibody (orb1262770)


Description	NEK3 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	This NEK3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 331-361 amino acids from the Central region of human NEK3.
Target	NEK3
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:2000
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	58 kDa
Uniprot ID	P51956
NCBI	P51956
Dilution Range	For WB starting dilution is: 1:2000
Expiration Date	12 months from date of receipt.



Western Blot
at 1:2000
dilution Lane
1: ...



Western Blot
at 1:2000
dilution Lane
1: ...



Western blot
analysis of
NEK3 using
rabb...