

Product Datasheet

PIP5K1C Antibody (orb1262791)

Description

PIP5K1C Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

This PIP5KI gamma (PIP5K1G) antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 637-668 amino acids from the C-terminal region of human PIP5KI gamma (PIP5K1G).

Target

PIP5K1C

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:50~100

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

73 kDa

Uniprot ID

[O60331](#)

NCBI


[O60331](#)

Dilution Range


For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:50~100

Expiration Date


12 months from date of receipt.



Western blot analysis of hPIP5K1G-L652 i...



Western blot analysis of PIP5K1C using P...



Formalin-fixed and paraffin-embedded hum...