

Product Datasheet

PIK3R5 Antibody (orb1262798)

Description

PIK3R5 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

This PI3KR5 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 392-421 amino acids from the Central region of human PI3KR5.

Target

PIK3R5

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

97 kDa

Uniprot ID

[Q8WYR1](#)

NCBI

[Q8WYR1](#)

Dilution Range

For WB starting dilution is: 1:1000

Expiration Date

12 months from date of receipt.

130

55
36
28
17

Western blot analysis of PI3KR5 Antibody...

72
55
43
34
26

(-) (+)

Western blot analysis of PIK3R5 using ra...