

Product Datasheet

PIK3CD Antibody (orb1262804)

Description

PIK3CD Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IF, IHC-P, WB

Applications

Immunogen

This PIK3CD antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 687-717 amino acids from the C-terminal region of human PIK3CD.

Target

PIK3CD

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For IF starting dilution is: 1:10~50For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

119 kDa

Uniprot ID

[O00329](#)

NCBI



[O00329](#)

Dilution Range


For IF starting dilution is: 1:10~50For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100

Expiration Date


12 months from date of receipt.



Fluorescent image of HeLa cell stained w...



Western blot analysis of anti-PIK3CD Pab...



Formalin-fixed and paraffin-embedded hum...