

Product Datasheet

PIK3C2B Antibody (orb1262810)

Description

PIK3C2B Antibody



Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated


Tested

IHC-P

Applications

Immunogen

This PI3KC2B antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 120-150 amino acids from the N-terminal region of human PI3KC2B.



Formalin-fixed and paraffin-embedded hum...

Target

PIK3C2B

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For IHC-P starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

185 kDa

Uniprot ID

[O00750](#)

NCBI

[O00750](#)

Dilution Range

For IHC-P starting dilution is: 1:10~50

Expiration Date

12 months from date of receipt.