

## Product Datasheet

### CYP3A43 Antibody (orb1262821)

## Description

CYP3A43 Antibody

### Species/Host

Rabbit

### Reactivity

Human

### Conjugation

Unconjugated

### Tested Applications

FC, IF, WB

### Immunogen

This CYP3A43 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 242-270 amino acids from the Central region of human CYP3A43.

### Target

CYP3A43

### Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

### Note

For research use only

### Application notes

For WB starting dilution is: 1:2000For FACS starting dilution is: 1:10~50For IF starting dilution is: 1:10~50

### Isotype

Rabbit Ig

### Clonality

Polyclonal

### MW

58 kDa

### Uniprot ID

[Q9HB55](#)

### NCBI


[Q9HB55](#)

### Dilution Range


For WB starting dilution is: 1:2000For FACS starting dilution is: 1:10~50For IF starting dilution is: 1:10~50

### Expiration Date


12 months from date of receipt.



Western Blot at  
1:2000 dilution +  
HepG2 ...



Flow cytometry  
analysis of 293  
cells (bo...



Confocal  
immunofluorescent  
analysis of C...