


Product Datasheet

PTGIS Antibody (orb1262822)

Description	PTGIS Antibody	72
Species/Host	Rabbit	43
Reactivity	Human	34
Conjugation	Unconjugated	26
Tested Applications	IHC-P, WB	
Immunogen	This CYP8A1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 472-500 amino acids from the C-terminal region of human CYP8A1.	Western blot analysis of CYP8A1 Antibody...
Target	PTGIS	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	57 kDa	
Uniprot ID	Q16647	
NCBI	Q16647	
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100	
Expiration Date	12 months from date of receipt.	



Formalin-fixed and paraffin-embedded hum...