

Product Datasheet

MYLK3 Antibody (orb1262836)

Description

MYLK3 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

This MYLK3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 40-69 amino acids from the N-terminal region of human MYLK3.

Target

MYLK3

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:50~100

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

88 kDa

Uniprot ID

[Q32MK0](#)

NCBI


[Q32MK0](#)

Dilution Range


For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:50~100

Expiration Date


12 months from date of receipt.



Western blot analysis of anti-MYLK3 Antibody



Western blot analysis of anti-MYLK3 Antibody



Formalin-fixed and paraffin-embedded human tissue