

Product Datasheet

PIM2 Antibody (orb1262853)

Description

PIM2 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, IP, WB

Applications

Immunogen

This PIM2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 277-308 amino acids from the C-terminal region of human PIM2.

Target

PIM2

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IP starting dilution is: 1:50~100For IHC-P starting dilution is: 1:50~100

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

34 kDa

Uniprot ID

[Q9P1W9](#)

NCBI


[Q9P1W9](#)

Dilution Range


For WB starting dilution is: 1:1000For IP starting dilution is: 1:50~100For IHC-P starting dilution is: 1:50~100

Expiration Date


12 months from date of receipt.



Western Blot
at 1:2000
dilution Lane
1: ...



Western blot
analysis of
anti-PIM2
Antib...



Western blot
analysis of
PIM2 using
rabb...