

Product Datasheet

CYP4F12 Antibody (orb1262872)

Description

CYP4F12 Antibody

NCI-H460
130

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

This CYP4F12 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 473-503 amino acids from the C-terminal region of human CYP4F12.

Target

CYP4F12

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

60 kDa

Uniprot ID

[Q9HCS2](#)

NCBI


[Q9HCS2](#)

Dilution Range

For WB starting dilution is: 1:1000

Expiration Date

12 months from date of receipt.



Western blot analysis of anti-CYP4F12 An...