

Product Datasheet

ALDH4A1 Antibody (orb1262886)

Description

ALDH4A1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

This ALDH4A1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 288-314 amino acids from the Central region of human ALDH4A1.

Target

ALDH4A1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For IHC-P starting dilution is: 1:50~100For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

62 kDa

Uniprot ID

P30038

NCBI


P30038

Dilution Range


For IHC-P starting dilution is: 1:50~100For WB starting dilution is: 1:1000

Expiration Date


12 months from date of receipt.



Formalin-fixed and paraffin-embedded hum...



Western blot analysis of anti-ALDH4A1 An...



Western blot analysis of ALDH4A1 using r...