

## Product Datasheet

### AMPD3 Antibody (orb1262891)

## Description

## AMPD3 Antibody

### Species/Host

Rabbit

### Reactivity

Human, Mouse

### Conjugation

Unconjugated

### Tested

WB

### Applications

### Immunogen

This AMPD3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 325-356 amino acids from the Central region of human AMPD3.

### Target

AMPD3

### Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

### Note

For research use only

### Application notes

For WB starting dilution is: 1:1000

### Isotype

Rabbit Ig

### Clonality

Polyclonal

### MW

89 kDa

### Uniprot ID

[Q01432](#)

### NCBI


[Q01432](#)

### Dilution Range


For WB starting dilution is: 1:1000

### Expiration Date


12 months from date of receipt.



Western blot analysis in human placenta ...



Western blot analysis in human placenta ...



Western blot analysis of anti-AMPD3 Anti...