

Product Datasheet

ALDH1A3 Antibody (orb1262906)

Description

ALDH1A3 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

This ALDH1A3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 318-348 amino acids from the Central region of human ALDH1A3.

Target

ALDH1A3

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

56 kDa

Uniprot ID

[P47895](#)

NCBI


[P47895](#)

Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western blot analysis of anti-ALDH1A3 Antibody



Western blot analysis of ALDH1A3 Antibody



Formalin-fixed and paraffin-embedded human tissue