

Product Datasheet

YARS2 Antibody (orb1262917)

Description

YARS2 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

This YARS2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 43-71 amino acids from the N-terminal region of human YARS2.

Target

YARS2

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

53 kDa

Uniprot ID

[Q9Y2Z4](#)

NCBI

[Q9Y2Z4](#)


Dilution Range

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50


Expiration Date

12 months from date of receipt.

K562
130



Western blot analysis of anti-YARS2 Anti...



Formalin-fixed and paraffin-embedded hum...