

Product Datasheet

ACVRL1 Antibody (orb1262933)

Description

ACVRL1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications

Immunogen

This ACVRL1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 38-68 amino acids from the N-terminal region of human ACVRL1.

Target

ACVRL1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50For IHC-P starting dilution is: 1:50~100

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

56 kDa

Uniprot ID

[P37023](#)

NCBI


[P37023](#)

Dilution Range


For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50For IHC-P starting dilution is: 1:50~100

Expiration Date


12 months from date of receipt.



Western blot of ACVRL1
Pab. TOP
LEFT: Mo...



Western blot analysis of ACVRL1 using ra...



Flow cytometric analysis of HepG2 cells

...