


Product Datasheet

ACE Antibody (orb1262938)

Description		ACE Antibody	 Western blot analysis of anti-ACE Antibody...
Species/Host		Rabbit	
Reactivity		Human, Mouse	
Conjugation		Unconjugated	
Tested Applications		WB	
Immunogen		This ACE antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 1274-1306 amino acids from the C-terminal region of human ACE.	
Target		ACE	
Preservatives		Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance		Liquid	
Concentration		batch dependent	
Storage		Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note		For research use only	
Application notes		For WB starting dilution is: 1:1000	
Isotype		Rabbit Ig	
Clonality		Polyclonal	
MW		150 kDa	
Uniprot ID		P12821	
NCBI		P12821	
Dilution Range		For WB starting dilution is: 1:1000	
Expiration Date		12 months from date of receipt.	