


Product Datasheet

HRAS Antibody (orb1262955)


Description

HRAS Antibody


Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	This HRAS antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 104-128 amino acids from the Central region of human HRAS.
Target	HRAS
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:500For IHC-P starting dilution is: 1:10~50For IF starting dilution is: 1:10~50
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	21 kDa
Uniprot ID	P01112
NCBI	P01112
Dilution Range	For WB starting dilution is: 1:500For IHC-P starting dilution is: 1:10~50For IF starting dilution is: 1:10~50
Expiration Date	12 months from date of receipt.



Western blot analysis of lysate from 293...



Western blot analysis of anti-HRAS Antib...



Formalin-fixed and paraffin-embedded hum...