

Product Datasheet

CDK5R1 Antibody (orb1262965)

Description

CDK5R1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

This CDK5R1(p35) antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 272-301 amino acids from the C-terminal region of human CDK5R1(p35).

Target

CDK5R1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For IHC-P starting dilution is: 1:10~50 For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

34 kDa

Uniprot ID

[Q15078](#)

NCBI


[Q15078](#)

Dilution Range


For IHC-P starting dilution is: 1:10~50 For WB starting dilution is: 1:1000

Expiration Date

12 months from date of receipt.



Formalin-fixed and paraffin-embedded hum...



Western blot analysis in A2058 cell line...