




Product Datasheet

PUM2 Antibody (orb1262969)

Description	PUM2 Antibody	
Species/Host	Rabbit	 Western blot analysis of PUM2 using PUM2...
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	FC, IHC-P, WB	
Immunogen	This PUM2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 159-189 amino acids from human PUM2.	
Target	PUM2	 Formalin-fixed and paraffin-embedded hum...
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	 Flow cytometric analysis of HepG2 cells
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	114 kDa	
Uniprot ID	Q8TB72	
NCBI	Q8TB72	
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50	
Expiration Date	12 months from date of receipt.	