

Product Datasheet

TRIB1 Antibody (orb1262973)

Description
TRIB1 Antibody
Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

This TRIB1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 341-370 amino acids from the C-terminal region of human TRIB1.

Target

TRIB1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:2000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

41 kDa

Uniprot ID
Q96RU8
NCBI
Q96RU8
Dilution Range


For WB starting dilution is: 1:2000

Expiration Date


12 months from date of receipt.



Western Blot
at 1:2000
dilution +
human ...



Western Blot
at 1:2000
dilution +
human ...



Western blot
analysis of
anti-TRIB1
Pab ...