




Product Datasheet

DDR2 Antibody (orb1262996)

Description	DDR2 Antibody	
Species/Host	Rabbit	 Western blot analysis of TYRO10 using ra...
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	FC, IHC-P, WB	
Immunogen	This TYRO10 antibody is generated from rabbits immunized with a his tag recombinant protein of human TYRO10.	
Target	DDR2	 Western blot analysis of TYRO10 Antibody...
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	 Formalin-fixed and paraffin-embedded hum...
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	97 kDa	
Uniprot ID	Q16832	
NCBI	Q16832	
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50	
Expiration Date	12 months from date of receipt.	