

Product Datasheet

FGFR2 Antibody (orb1263052)

Description

FGFR2 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IF, IHC-P, WB

Applications

Immunogen

This FGFR2 antibody is generated from rabbits immunized with a his tag recombinant protein of human FGFR2.

Target

FGFR2

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For IF starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

92 kDa

Uniprot ID

P21802

NCBI


P21802

Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For IF starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western blot analysis of FGFR2 using rabbit antibody...



Western blot analysis of FGFR2 Antibody ...



Formalin-fixed and paraffin-embedded human tissue...