



Product Datasheet

FGFR1 Antibody (orb1263056)

Description	FGFR1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	This FGFR1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 132-161 amino acids from human FGFR1.
Target	FGFR1
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:1000
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	92 kDa
Uniprot ID	P11362
NCBI	P11362
Dilution Range	For WB starting dilution is: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of anti-FGFR1 Antibody



Western blot analysis of FGFR1 using rabbit anti-FGFR1 antibody