

Product Datasheet

POLDIP2 Antibody (orb1263077)

Description

POLDIP2 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

FC, WB

Applications
Immunogen

This POLDIP2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 320-349 amino acids from the C-terminal region of human POLDIP2.

Target

POLDIP2

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


42 kDa

Uniprot ID
[Q9Y2S7](#)
NCBI
[Q9Y2S7](#)
Dilution Range


For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western blot analysis of POLDIP2 antibody...



Flow cytometry analysis of MCF-7 cells (...)



Western blot analysis of POLDIP2 using r...