

Product Datasheet

EPHA1 Antibody (orb1263102)

Description

EPHA1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications
Immunogen

This EphA1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 29-60 amino acids from the N-terminal region of human EphA1.

Target

EPHA1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

108 kDa

Uniprot ID
P21709
NCBI
P21709
Dilution Range

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50


Expiration Date

12 months from date of receipt.


250

 75
50
37
25
20

Western blot analysis of anti-EphA1 N-te...



Formalin-fixed and paraffin-embedded hum...



Flow cytometric analysis of HepG2 cells

...