

Product Datasheet

YARS Antibody (orb1263125)

Description

YARS Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

This Tyrosyl tRNA synthetase (YARS) antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 450-479 amino acids from the C-terminal region of human Tyrosyl tRNA synthetase (YARS).

Target

YARS

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

59 kDa

Uniprot ID
P54577
NCBI
P54577
Dilution Range

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50


Expiration Date

12 months from date of receipt.

250
130

55
36
28
17

Western blot analysis of anti-YARS Antib...



Formalin-fixed and paraffin-embedded hum...